

# Set Up MastaCube

Connection: Ethernet 

## 1. Included in box



1 Power cable



1 Ethernet cable



1 MastaCube



## 2. Set up



**Connect Power & Ethernet**



**Optional:**  
Attach a USB keyboard  
(not included)



**Turn power ON**  
at wall and wait for  
MastaCube to initiate

## 3. Set Language



**Scroll** to your preferred language using the grey buttons.  
**Hold** (2s) the right button to select.  
**Press** the right button to confirm the change.

## 4. Register

Scan QR on MastaCube screen and complete registration

**OR**

In a web browser, navigate to [www.mastatest.com/register](http://www.mastatest.com/register) and complete registration



*Hint: The Serial Number is on the bottom of the MastaCube, and the Hardware ID is displayed on the MastaCube screen.*

## 5. Start testing

Scan for a detailed guide on how to start a Mastatest.



[www.mastatest.com/mc](http://www.mastatest.com/mc)

